

# High specification patient trolleys



## EMED 1600 & EMED 1602 PATIENT TROLLEYS




E-MED/1600 (Fixed height trolley with independent brake)  
E-MED/1602 (Fixed height trolley with constant brake)



# E-MED 1600 & 1602 Fixed height trolleys

**E-MED 1602**  
(Constant brake)




Brake is constantly applied - Push down on bar to release

### Key features

- The E-Med 1600 & 1602 fixed height patient trolleys feature two separate platform sections:
- Independent backrest section (angle 0-80°)
- Femoral & leg rest section (fixed position)
- Fixed height 2 section patient trolley
- The trolley has a fixed height of 78cm
- Backrest with multiple positions for patient comfort can be positioned up to 80°
- Lightweight frame with side bumpers enables ease of trolley movement
- 200mm ø castors as standard
- A range of optional push handles and push bars available
- Compass and folding side rails available
- Variety of accessories can be used with the trolley
- Compliant to 93/42/EEC Class 1 Medical Device

**E-MED 1600**  
(Independent brake)



FIXED HEIGHT OF 78CM

Product Dimensions	
Total length	204cm
Total width	70cm
Platform surface area	198 x 58cm
Trolley height	78mm
Backrest angle	80°
Mattress	193 x 58 x 5cm
Total weight	48kg
Max. patient weight	265kg (42st)
Safe working load	300kg (47st)

Product:	Code:	Description:
E-Med 1600	EMED/1600	Fixed height - independent brake
E-Med 1602	EMED/1602	Fixed height - constant brake

The E-MED 1600 & 1602 come complete with:  
4 Castors 200mm ø / 2 section trolley (backrest, pelvic/femoral/legrest section / Side bumpers

## Accessories

Optional accessories:	
EMED/ACC/Z1501	Antistatic mattress for fixed height - constant brake (1602 only)
EMED/ACC/Z1506	Antistatic mattress for fixed height - independent brake (1600 only)
EMED/ACC/Z1524	Adjustable IV pole with 4 hooks
EMED/ACC/Z1525	Paper roll holder
EMED/ACC/Z1520	Oxygen bottle holder
EMED/ACC/Z1609	Clothes basket
EMED/ACC/S0070	Drainage bag holder (not shown below)

Factory fitted accessories:	
EMED/ACC/Z1604	Folding side rails
EMED/ACC/Z1607	Compass side rails (not shown below)
EMED/ACC/Z1565	Radio transparent surface (not shown below)
EMED/ACC/Z1535	Push handle
EMED/ACC/Z1518	Monitor support (not shown below)

## Trolley kit packages

Fixed height trolley: Independent brake KIT CODE: EMED/1600/KIT	
EMED/1600	E-Med 1600 Trolley
EMED/ACC/Z1604	Folding side rails
EMED/ACC/Z1506	Antistatic mattress for fixed height independent brake
EMED/ACC/Z1524	IV stand with 4 hooks
EMED/ACC/Z1520	Oxygen bottle holder
EMED/ACC/Z1525	Paper roll holder
EMED/ACC/Z1609	Clothes basket

Fixed height trolley: Constant brake KIT CODE: EMED/1602/KIT	
EMED/1602	E-Med 1602 Trolley
EMED/ACC/Z1604	Folding side rails
EMED/ACC/Z1501	Antistatic mattress for fixed height constant brake
EMED/ACC/Z1524	IV stand with 4 hooks
EMED/ACC/Z1520	Oxygen bottle holder
EMED/ACC/Z1525	Paper roll holder
EMED/ACC/Z1609	Clothes basket


## Service centres nationwide

With service centres strategically located nationwide, we guarantee an optimum response, supported by our team of field service engineers. Their experience and depth of product knowledge ensures equipment is put back into service promptly, efficiently and cost effectively.

We carry comprehensive stocks of spare parts both for the Sidhil range and for other manufacturers' products at our service centres. We also provide decontamination services to protocols as advised by a clinical scientist in infection control, staffed by trained personnel operating in accordance with MHRA DB2006 (05).

Our service centre processes are embedded within Sidhil's established quality standards to ISO 9001 levels.

For more information on our service and maintenance capabilities, contact the Service Team on: 01422 233 136.



- Glasgow
- Sunderland
- Halifax
- Warrington
- Kettering
- Orpington
- N. Ireland

## Service & Maintenance & Rentals

Product Servicing - All Makes / Decontamination / Spares  
Short term solutions available as an alternative to outright purchase



Sidhil Business Park  
Holmfild, Halifax,  
West Yorkshire  
HX2 9TN

Tel: +44 (0) 1422 233000  
Fax: +44 (0) 1422 233010  
Email: [sales@sidhil.com](mailto:sales@sidhil.com)  
Web: [www.sidhil.com](http://www.sidhil.com)



Independently verified to IEC 60601-1 and IEC 60601-2-52 through a 3rd party UKAS approved test house