

Hebden nursing and residential care bed range...



Hebden Nursing Care Bed

Further expanding the company's range of electrically operated, fully profiling beds developed specifically for nursing and residential care home environments, healthcare specialists Sidhil have introduced five new scissor action nursing care beds designed around the developing requirements of the care home market.

The all-new Hebden beds provide a comprehensive portfolio of functions developed specifically to meet the needs of nursing and residential care home applications, complemented by an attractive domestic appearance

with high quality covered wooden head and foot ends available in a range of wood colour options to match interior décor schemes.




Key features

- Electrically operated, variable height fully profiling nursing care bed
- Variable height mattress platform 27.5 - 85cm for ease of client access, reducing the need for lifting and handling
- Full length integral wooden side rails can be released at both head and foot ends for safety and convenience
- Underbed clearance of 29cm to facilitate hoist access is not impeded
- Infinitely variable tilt function: head down tilt from 0 - 15°, foot down tilt from 0 - 18°
- Lockable handset
- All bed sections weigh less than 25kg in line with Health & Safety specifications
- Easy to move, assemble and disassemble even where space for manoeuvre is limited
- Maximum user weight 180kg (28st)
- 3 year warranty (1st year parts & labour, 2nd & 3rd years parts only)
- 10 year warranty available with Service & Maintenance Contract

Customise your Hebden bed

Wood finish options



Choose from 10 classic wood finishes. Standard finish is Light Oak

*Please specify wood colour when ordering

Hebden Accessories



Cantilever bed cradle



Adjustable backrest



Hebden side rail pads



Hebden side rail mesh pads



Lifting pole



Grab handle

Accessories	
Lifting Pole	1275/LP
Cantilever Bed Cradle	1369
9 Button Handset	CONTROL/121
Side Rail Channel Kit <small>(for use when the bed does not have side rails)</small>	1275/CHANNEL/ INFIL/KIT
Sides	
Hebden Side Rail Pads	1316/PS/74/CREAM
Hebden Side Rail Mesh Pads	1275/PS
Grab Handle	
Grab Handle	1214/GRAB

Hebden nursing and residential care bed range...

Bed Dimensions

Height 27.5 to 85 (to mattress platform)	Width 105cm (Incl. side rails)
27.5 to 85 (to mattress platform) Length	Length 208cm
Maximum User Weight	190kg / 28st

Product codes

Hebden Bed 01 (without side rails):	HEB/01/125/ LOAK
Hebden Bed 01 (with side rails):	HEB/01/ SR/125/LOAK
Hebden Bed 02 (without side rails):	HEB/02/125/ LOAK
Hebden Bed 02 (with side rails):	HEB/02/ SR/125/LOAK
Hebden Bed 03 (without side rails):	HEB/03/125/ LOAK
Hebden Bed 03 (with side rails):	HEB/03/ SR/125/LOAK
Hebden Bed 04 (without side rails):	HEB/04/125/ LOAK
Hebden Bed 04 (with side rails):	HEB/04/ SR/125/LOAK
Hebden Bed 05 (without side rails):	HEB/05/125/ LOAK
Hebden Bed 05 (with side rails):	HEB/05/ SR/125/LOAK




Service centres nationwide

With service centres strategically located nationwide, we guarantee an optimum response, supported by our team of field service engineers. Their experience and depth of product knowledge ensures equipment is put back into service promptly, efficiently and cost effectively.

We carry comprehensive stocks of spare parts both for the Sidhil range and for other manufacturers' products at our service centres. We also provide decontamination services to protocols as advised by a clinical scientist in infection control, staffed by trained personnel operating in accordance with MHRA DB2006 (05).

Our service centre processes are embedded within Sidhil's established quality standards to ISO 9001 levels.

For more information on our service and maintenance capabilities, contact the Service Team on: 01422 233 136.



- Glasgow
- Sunderland
- Halifax
- Warrington
- Kettering
- Orpington
- N. Ireland

Service & Maintenance & Rentals

Product Servicing - All Makes / Decontamination / Spares
Short term solutions available as an alternative to outright purchase



Sidhil Business Park
Holmfield, Halifax,
West Yorkshire
HX2 9TN

Tel: +44 (0) 1422 233000
Fax: +44 (0) 1422 233010
Email: sales@sidhil.com
Web: www.sidhil.com



Independently verified to IEC 60601-1 and IEC 60601-2-52
through a 3rd party UKAS approved test house