

# HYBRID-AIR PLUS

## Mattress

“A non-powered hybrid suitable for the heavier individual”



### Features Include:

- › Castellated side formers to prolong longevity
- › Visco-elastic foam covers the mattress surface – enhancing comfort
- › White inner substrate – fluid ingress becomes easily identifiable to aid infection control guidelines
- › Non-slip base – prevents the mattress from sliding out of position, giving reassurance to the carer
- › Sodium Hypochlorite resistant cover – withstands a higher concentration of cleaning substances
- › 360° zip to prolong longevity with high frequency welded seams
- › Fully welded product to protect against fluid ingress

### Specifications:

CODES: 1003300002

Mattress Height (Inflated):  
2000mm

Mattress Width (Inflated):  
1200mm

Mattress Depth (Inflated):  
200mm

Mattress Warranty:  
2 years

Fire Retardancy:  
Crib 5

Risk Factor:  
Very High Risk

Max Weight:  
319kg / 50st

The Hybrid-Air Plus is a support surface that combines high specification foam and air cells which adjust with body movement. The support surface caters for the heavier individual and sits within a 'U' shaped foam base, providing maximum support.

Leaflet Code: LL326

Tel: 0845 0600 333 / + (0)1274 475000 | Fax: 0845 0600 334

Email: [info@drivedevilbiss.co.uk](mailto:info@drivedevilbiss.co.uk) | Web: [www.drivedevilbiss.co.uk](http://www.drivedevilbiss.co.uk)

**drive**  
DeVilbiss<sup>®</sup>  
HEALTHCARE