

# HYBRID-POWER

## Cushion

“Visco-elastic foam combined with alternating cells”

### Features Include:

- › Compact and light weight power unit
- › Cover can withstand a higher concentration of cleaning substances and offers greater resistance to abrasion
- › High frequency welded seams protect against fluid ingress
- › White inner substrate - fluid ingress becomes easily identifiable to aid infection control guidelines
- › Automatic lock out feature
- › Low pressure alarm and indicator
- › Weight limit of 30 stones

CODE: 1002290001



The Hybrid-Power Cushion combines the pressure redistribution properties of high specification foam, with alternating pressure relief. The visco-elastic foam conforms to the shape of the body, while the alternating cells assist with pressure relief when the patients risk level changes.

Leaflet Code: LL251

[www.drivedevilbiss.co.uk](http://www.drivedevilbiss.co.uk)

**drive**  
**DeVilbiss**  
HEALTHCARE