

Solite Pro
Introducing Sidhil's latest
generation of community beds...

SOLITE  **PRO**



drive
 **DeVilbiss**
HEALTHCARE

SIDHIL

Solite Pro community beds

Acknowledged to be the highest specification product in its class, Sidhil's Solite range has been constantly developed and improved over the years to provide the optimum combination of client comfort, versatility, functionality and performance.

The new generation of Solite Pro community beds brings to market high quality, UK manufactured electrically profiling beds with electric backrest and knee break function and variable height control operated via improved actuators which add to the strength of the bed. New integral full length side rails conforming to all gapping

requirements feature a unique central 'stepped down' profile for easier ingress and egress.

Solite Pro beds come complete with a five year warranty as standard, which can be extended to ten years for applications where a Sidhil bed service and maintenance contract remains in place.



Verified to IEC 60601-1 General Medical Electrical Equipment standard, IEC 60601-2-52 Medical Bed standard and IEC 60601-1-11 Medical Electrical Equipment standard specific to homecare requirements, the Solite Pro is a complete redesign of the Solite UK, now offering additional functionality whilst retaining an attractive domestic look and feel.




Solite Pro

Key Features

- Redesigned platform for improved patient comfort when profiling
- Individual sections strengthened for extended performance and reliability
- Optional integral full length side rails for safety and ease of use
- Improved integral leg extension with safety features to prevent entrapment
- Improved auto regressing backrest prevents pinching in the sacral area
- Infinitely variable foot down tilt - head down tilt available as an optional extra. New handset (CONTROL/038) would be required
- Four braked castors for manoeuvrability and safety
- Easy-use handset controls all functions - includes lockout function
- All individual sections of the bed are below 25kg, aiding assembly
- 5 year warranty as standard (can be extended to 10 years where a Sidhil service contract remains in place)
- Supplied complete with transport stand
- Maximum user weight 191kg (30 stone)
- Maximum safe working load 226kg (35.5 stone)

Bed ordering process



Bed Dimensions

	Height (cm)	Width (cm)	Length (cm)	Weight (kg)
Solite Pro	39.5 - 81.5 (to top of mattress platform)	96.5 (101.5 with side rail channels)	223 (239.5 when fully extended)	85
Solite Pro Low	29.5 - 71.5 (to top of mattress platform)	96.5 (101.5 with side rail channels)	223 (239.5 when fully extended)	85

Solite Pro introducing Sidhil's latest generation of community beds...

Beds	
Solite Pro Bed	SOLITE/PRO
Solite Pro Low Bed	SOLITE/PRO/LOW
Solite Pro Low with Battery Back Up	SOLITE/PRO/LOW/BB

Side Rail Channels	SOLITE/PRO/ADD/CHANNELS
Side Rail Channels - Factory fitted to bed	SOLITE/PRO/CHANNEL/FITTING

Side Rails	
Solite Pro Full Length Side Rails*	SOLITE/PRO/SR/A2
Solite Pro Extended Length Side Rails*	SOLITE/PRO/SR/EXT
Solite Pro Side Rail Height Extensions	SOLITE/PRO/SR/HE
Grange Side Rails	GRANGE/SR

* Please note: To fit this accessory to the Solite Pro then side rail channels must be present. Accessories should be fitted by a competent technician

Side Rail Pads and End Pads	
Solite Pro Standard Side Rail Pads	SOLITE/PRO/SR/PAD
Solite Pro Side Rail Length Extension Pads	SOLITE/PRO/SR/EXT/PAD
Solite Pro Side Rail Height Extension Pads	BED01/SR/HEXT/PAD
Solite Pro Standard Side Rail Mesh Cradle	SOLITE/PRO/EPS
General Side Rail Pads (for use with Grange Side Rails)	1316/PS/CREAM
Solite Pro End Pads	SOLITE/PRO/HF/PAD

Accessories	
Transfusion Pole	SOLITE/DP
Lifting Pole	SOLITE/LASER/LP/1
Cromer Cantilever Bed Cradle	1369
Spare Lockable Handset	CONTROL/127
Lockable Handset with Head Down Tilt Function	CONTROL/038
Solite Pro Side Rail Transport Accessory	SOLITE/PRO/SR/TS
Handset Locking Key	CONTROL/049

Grab Handles	
Elland Standard Grab Handle	1210/GRAB/MAG
Elland Short Grab Handle	1216/GRAB

Service centres nationwide

With service centres strategically located nationwide, we guarantee an optimum response, supported by our team of field service engineers. Their experience and depth of product knowledge ensures equipment is put back into service promptly, efficiently and cost effectively.


We carry comprehensive stocks of spare parts both for the Sidhil range and for other manufacturers' products at our service centres. We also provide decontamination services to protocols as advised by a clinical scientist in infection control, staffed by trained personnel operating in accordance with MHRA DB2006 (05).

Our service centre processes are embedded within Sidhil's established quality standards to ISO 9001 levels.

For more information on our service and maintenance capabilities, contact the Service Team on: 01422 233 136.

Service & Maintenance & Rentals

Product Servicing - All Makes / Decontamination / Spares
Short term solutions available as an alternative to outright purchase



- Glasgow
- Sunderland
- Halifax
- Warrington
- Kettering
- Orpington
- N. Ireland



Sidhil Business Park
Holmfild, Halifax,
West Yorkshire
HX2 9TN

Tel: +44 (0) 1422 233000
Fax: +44 (0) 1422 233010
Email: sales@sidhil.com
Web: www.sidhil.com



93/42/EEC

Certificate No. FM14550

Independently verified to IEC 60601-1 and IEC 60601-2-52 through a 3rd party UKAS approved test house